

FACS Analysis Checklist

Please follow the parameters below for bringing materials for flow cytometry.

- Two million cells per sample before labeling.
- One sample unstained cells.
- One control sample per stain.
- 5 ml Falcon 12 x 75 polystyrene tubes (BD Labware Cat #352054).
- Final volume should be 250 μ l to 500 μ l at approximately 2 million cells per sample.
- If your schedule will delay your arrival for more than 15 minutes, please give us a call.

Thank you for your attention to this information. We appreciate your business and look forward to working with you.